



**Haley's Joy®**  
Swing Therapy Anywhere!

## Mat Padding Instruction Manual

### WARNING

**READ ENTIRE MANUAL BEFORE USE. THIS SWING IS NOT A TOY.**

**THIS SWING IS ONLY TO BE USED BY TRAINED PERSONNEL, SUCH AS AN OCCUPATIONAL THERAPIST, PHYSICAL THERAPIST, SPEECH THERAPIST, ETC.**

**THIS SWING IS DESIGNED AS A THERAPEUTIC TOOL FOR A CHILD WITH DEVELOPMENTAL NEEDS. THE THERAPIST IS TO SUPERVISE AT ALL TIMES. USE OF THE SWING SHOULD ALWAYS BE IN A CONTROLLED FASHION WITH THE TRAINED PERSONNEL'S HAND ON THE SWING, GUIDING THE MOTION.**

**THE SWING SHOULD ALWAYS BE GUIDED BY PLACING THE HANDS ON THE SWING, ITSELF. PUSHING ON THE CHAINS COULD RESULT IN EARLY FAILURE OF THE CHAINS AND VINYL.**

**ONLY UTILIZE THE SWING WITH ONE INDIVIDUAL AT A TIME, AND ALL HANDS AND FEET SHOULD BE KEPT OUT FROM UNDER THE SEAT.**

**THIS SWING IS A LOW MOTION DEVICE; KEEP THE SWING'S OCCUPANT / SEAT WITHIN THE PERIMETER OF THE FOUR LEGS. FAILURE TO DO SO CAN RESULT IN TIPPING OF THE SWING.**

**NEVER ALLOW THE INDIVIDUAL TO STAND ON / IN THE SWING; EVEN DURING THERAPY.**

**THE SWING IS DESIGNED TO BE USED CLOSE TO THE GROUND.**

**THERE IS NO SEAT BELT. CARE MUST BE TAKEN NOT TO ALLOW THE CHILD TO FALL OFF THE SWING.**

**CHOKING HAZARD: DO NOT ALLOW EATING OR DRINKING ON THE SWING!**

**STRANGULATION HAZARD: DO NOT ATTACH STRINGS OR ROPES TO THE SWING OR PLACE THE SWING NEXT TO ANY STRINGS OR ROPES.**

**UNSUPERVISED USE OF THIS SWING CAN RESULT IN SEVERE INJURY. WHEN NOT IN USE, THE SWING MUST BE DISASSEMBLED AND STORED IN A SAFE PLACE, ONLY ACCESSIBLE TO TRAINED PERSONNEL.**



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please keep this Instruction Sheet in an accessible location. The Instruction Sheet must be read and familiar to anyone who is using the Mat Padding. If this Instruction Sheet is misplaced, please visit the product page on the FlagHouse website or email ([haleysjoy@flaghouse.com](mailto:haleysjoy@flaghouse.com)) for a replacement.

## IMPORTANT

Please check equipment immediately for any damage that may have occurred in shipping. Contact FlagHouse immediately if there is a damaged piece 800.793.7900 (US & International) or 800.265.6900 (Canada).

## MAT PADDING ASSEMBLY INSTRUCTIONS


Mat padding is provided, as an option, for each size swing. The number of components that are supplied for each swing and how they are to be laid out are shown below. Just align the 'teeth' and press them together. Although either side of the pad could go up or down, the smooth side probably provides the best grip on hard surfaces. Align the 'teeth' of the edges pieces and press them onto the pads.

The longer direction of the On the Go - Size 2 is aligned with the long dimensions of the On the Go - Size 2 Swing (in line with the backbone).

Note: Align the On the Go Swings to be in the center of the mat at all times.

### The On the Go - Size 1 Swing:


The On the Go - Size 1 has 9 tiles, 4 corner pieces and 8 edge pieces. The layout is as follows:



approximately 70" x 70" without edges  
approximately 73" x 73" with edges

### The On the Go - Size 2 Swing:


The On the Go - Size 2 has 12 tiles with 4 corner edges and 10 edge pieces. The layout is as follows:



approximately 70" x 93 ½" without edges  
approximately 73" x 96 ½" with edges

### The On the Go - Size 3 Swing:

The On the Go - Size 3 has 16 tiles with 4 corner edges and 12 edge pieces. The layout is as follows:



approximately 93 ½" x 93 ½" without edges  
approximately 96 ½" x 96 ½" with edges



601 FlagHouse Drive | Hasbrouck Heights, NJ 07604 | 800.793.7900 | [flaghouse.com/haleysjoy](http://flaghouse.com/haleysjoy)



235 Yorkland Boulevard | North York, ON M2J 4Y8 | 800.265.6900 | [flaghouse.ca/haleysjoy](http://flaghouse.ca/haleysjoy)

